

Medline Sterile Procedure Nitrile

Sterile nitrile - powder-free examination gloves



Reference(s)

MDS2294 (S: 6/6.5)
MDS2295 (M: 7/7.5)
MDS2296 (L: 8/8.5)
MDS2297 (XL: 9)

Specifications

Material: Nitrile

Cuff: Beaded

Latex free: Yes

Accelerator: Carbamate

Grip:



Colour:



Thickness

Fingertip: 0.12 mm

Palm: 0.10 mm

Cuff: 0.06 mm

Dimensions

Box: 150 mm W x 250 mm L x 135 mm H

Case: 280 mm W x 460 mm L x 260 mm H

Packaging

50 gloves per box,

6 boxes per case



Technical sheet

Medline Sterile Procedure Nitrile

| | |
|-----------------------------------|--|
| Primary material | Nitrile |
| Colour | Blue |
| Grip | Textured fingertips |
| Cuff design | Beaded |
| Latex free | Yes |
| Accelerator | Carbamate |
| Thickness | Finger tip 0.12 mm Palm 0.10 mm Cuff 0.06 mm |
| Dimensions | Box: 150 mm W x 250 mm L x 135 mm H Case: 280 mm W x 460 mm L x 260 mm H |
| Packaging | 50 gloves per box 6 boxes per case |
| Physical properties | <p>EN 455-1 Freedom From Holes: 1.0 AQL (standard requests compliance with AQL 1.5)</p> <p>EN 455-2 Dimensions: Length: 302 mm Width: 106 mm</p> <p>Force at Break: Before challenge: 7.8 N After challenge: 8.3 N</p> <p>EN 455-3 Powder Content: 0.02 milligrams per glove</p> <p>EN 455-4 Shelf Life Properties: Meets requirements and testing for shelf-life determination of 35 months</p> |
| Chemical resistance | The resistance to chemicals has been assessed in accordance to EN 16523-1 and chemo drugs under ASTM D 6978. Results can be obtained upon request. |
| Viral penetration | Tested and passed, in accordance with EN374-5. |
| Regulations and quality standards | <p>Class 1 non-sterile Medical Device per Directive 93/42/EEC.</p> <p>Complies with standards EN 455-1, 2, 3 and 4.</p> <p>Medline manufacturing locations are certified to EN ISO 13485.</p> |
| PPE Certification | Under the requirements of Personal Protective Equipment Regulation (EU)2016/425 Category III. Complies with standards EN 420, EN ISO 374-1,2,4,5 and EN 16523-1. |
| Storage recommendations | Protect from freezing. Avoid excessive heat. Keep dry. Product should be shielded from direct sunlight, fluorescent lighting, X-rays, moisture and ozone. |



For more information, contact your Medline Account Manager or call 1800 110 511.

**Together,
advancing
healthcare**



Medline International Two Australia Pty Ltd
2 Fairview Place, Marsden Park NSW 2765
1800 110 511
medline.com.au