

MediGuard® Special Nitrile

Non-sterile nitrile - powder-free - examination gloves



Reference(s)

MGSNXS (XS)
MGSNS (S)
MGSNM (M)
MGSNL (L)
MGSNXL (XL)


Specifications

Material: Nitrile

Cuff: Beaded

Latex free: Yes

Accelerator: Carbamate



Thickness

Fingertip: 0.10 (\pm 0.03) mm

Palm: 0.08 (\pm 0.03) mm

Cuff: 0.07 (\pm 0.03) mm

Dimensions

Box: 125 mm W x 245 mm L x 74 mm H

Case: 257 mm W x 382 mm L x 255 mm H

Packaging

200 gloves per box, (XL - 180 gloves per box)

10 boxes per case



Size

X-small

Yellow

Small

Orange

Medium

Blue

Large

Green

X-large

Black

Technical sheet

MediGuard® Special Nitrile

Primary material	Nitrile
Colour	Violet blue
Grip	Textured fingertips
Cuff design	Beaded
Latex free	Yes
Accelerator	Carbamate
Thickness	Finger tip 0.10 (± 0.03) mm Palm 0.08 (± 0.03) mm Cuff 0.07 (± 0.03) mm
Dimensions	Box: 125 mm W x 245 mm L x 74 mm H Case: 257 mm W x 382 mm L x 255 mm H
Packaging	200 gloves per box, (XL - 180 gloves per box) 10 boxes per case
Physical properties	EN 455-1 Freedom From Holes: 1.0 AQL (standard requests compliance with AQL 1.5) EN 455-2 Dimensions: Length: 240 mm Width: 97 mm Force at Break: Before challenge: 7.0 N After challenge: 7.5 N EN 455-3 Powder Content: 0.20 milligrams per glove EN 455-4 Shelf Life Properties: Meets requirements and testing for shelf-life determination of 35 months
Chemical resistance	The resistance to chemicals has been assessed in accordance to EN 16523-1 and chemo drugs under ASTM D 6978.
Viral penetration	Tested and passed, in accordance with EN374-5
Regulations and quality standards	Class 1 non-sterile Medical Device per Regulation (EU) 2017/745, and UK Medical Device Regulation. Complies with standards EN 455-1, 2, 3 and 4. Medline manufacturing locations are certified to EN ISO 13485
Food Contact	Tested according to European Regulation (EU) No 10/2011
PPE Certification	Under the requirements of Personal Protective Equipment Regulation (EU)2016/425 Category III. Complies with standards EN 420, EN ISO 374-1,2,4,5 and EN 16523-1.

Storage recommendations

Protect from freezing. Avoid excessive heat. Keep dry. Product should be shielded from direct sunlight, fluorescent lighting, X-rays, moisture and ozone.



For more information, contact your Medline Account Manager or call 1800 110 511.

**Together,
advancing
healthcare**



Medline International Two Australia Pty Ltd
2 Fairview Place, Marsden Park NSW 2765
1800 110 511
medline.com.au