

SensiCare® Free

Non-sterile powder-free and accelerator-free nitrile examination gloves




Reference(s)

MDSFREE-XS (XS)
MDSFREE-S (S)
MDSFREE-M (M)
MDSFREE-L (L)
MDSFREE-XL (XL)

Specifications

Material: Nitrile
Cuff: Beaded
Latex free: Yes
Accelerator: None



Thickness

Fingertip: 0.09 mm
Palm: 0.07 mm
Cuff: 0.05 mm

Dimensions

Box: 240 mm W x 120mm L x 77 mm H
Case: 400 mm W x 248 mm L x 255 mm H

Packaging

200 gloves per box, (XL - 180 gloves per box)
10 boxes per case



Size

X-small	Small	Medium	Large	X-large
Yellow	Orange	Blue	Green	Black

Technical sheet

SensiCare® Free

Primary material	Nitrile
Colour	Blue
Grip	Finger-tip textured
Cuff design	Beaded
Latex free	Yes
Accelerator	None
Thickness	Finger tip 0.09 mm Palm: 0.07 mm Cuff 0.05 mm
Dimensions	Box: 240 mm W x 120 mm L x 77 mm H Case: 400 mm W x 248 mm L x 255 mm H
Packaging	200 gloves per box, (XL - 180 gloves per box) 10 boxes per case
Physical properties	EN 455-1 Freedom From Holes: 1.0 AQL EN 455-2 Dimensions: Length: 243 mm Width: 98 mm Force at Break: Before challenge: 8.2 N After challenge: 7.9 N EN 455-3 Powder Content: 0.10 milligrams per glove EN 455-4 Shelf Life Properties: Meets requirements and testing for shelf-life determination of 35 months
Chemical resistance	The resistance to chemicals has been assessed in accordance with EN 16523-1 and chemo drugs under ASTM D 6978. Results can be obtained upon request
Viral penetration	Tested and passed, in accordance with EN374-5
Regulations and quality standards	Class 1 non-sterile Medical Device per Regulation (EU) 2017/745. Complies with standards EN 455-1, 2, 3 and 4. Medline manufacturing locations are certified to EN ISO 13485
Food contact	Tested according to European Regulation (EU) No 10/2011
PPE Certification	Under the requirements of Personal Protective Equipment Regulation (EU)2016/425 Category III. Complies with standards EN 420, EN ISO 374-1,2,4,5 and EN 16523-1
Storage recommendations	Protect from freezing. Avoid excessive heat. Keep dry. Product should be shielded from direct sunlight, fluorescent lighting, X-rays, moisture and ozone.



For more information, contact your Medline Account Manager or call 1800 110 511.

**Together,
advancing
healthcare**



Medline International Two Australia Pty Ltd
2 Fairview Place, Marsden Park NSW 2765
1800 110 511
medline.com.au