

Aloetouch® Extended Cuff

Non-sterile nitrile - powder-free - ambidextrous examination gloves with aloe vera



Reference(s)

MDS195184 (S)
MDS195185 (M)
MDS195186 (L)
MDS195187 (XL)

Specifications

Material: Nitrile

Cuff: Beaded

Latex free: Yes

Accelerator: Carbamate

Grip:



Colour:



Thickness

Fingertip: 0.08 mm

Palm: 0.07 mm

Cuff: 0.09 mm

Dimensions

Box: 120 mm W x 270 mm L x 58 mm H

Case: 255 mm W x 310 mm L x 290 mm H

Packaging

50 gloves per box

10 boxes per case



Size

X-small

Small

Medium

Large

X-large

Yellow

Orange

Blue

Green

Black

Technical sheet

Aloetouch® Extended Cuff

Primary material	Nitrile
Colour	Green
Grip	Full textured
Cuff design	Beaded
Latex free	Yes
Accelerator	Carbamate
Thickness	Finger tip 0.08 mm Palm 0.07 mm Cuff 0.09 mm
Dimensions	Box: 120 mm W x 270 mm L x 58 mm H Case: 255 mm W x 310 mm L x 290 mm H
Packaging	50 gloves per box, 10 boxes per case
Physical properties	EN 455-1 Freedom From Holes: 1.5 AQL EN 455-2 Dimensions: Length: 307 mm Width: 98 mm Force at Break: Before challenge: 8.6 N After challenge: 8.0 N EN 455-3 Powder Content: 0.02 milligrams per glove EN 455-4 Shelf Life Properties: Meets requirements and testing for shelf-life determination of 35 months
Chemical resistance	The resistance to chemicals has been assessed in accordance with EN 16523-1. Results can be obtained upon request
Viral penetration	Tested and passed, in accordance with ASTM F1671 standard
Regulations and quality standards	Class 1 non-sterile Medical Device per Regulation (EU) 2017/745 Complies with standards EN 455-1, 2, 3 and 4. Medline manufacturing locations are certified to EN ISO 13485
Storage recommendations	Protect from freezing. Avoid excessive heat. Keep dry. Product should be shielded from direct sunlight, fluorescent lighting, X-rays, moisture and ozone.



For more information, contact your Medline Account Manager or call 1800 110 511.

**Together,
advancing
healthcare**



Medline International Two Australia Pty Ltd
2 Fairview Place, Marsden Park NSW 2765
1800 110 511
medline.com.au