

Central Sterile Supply Department (CSSD) Products



**Your partner
in sterilisation
and beyond**

About Medline

We're a global healthcare company dedicated to being easy to partner with, solving problems, and providing good value. In Australia, from our network of warehouses nationally, we manufacture and distribute an expansive portfolio of products for hospitals, healthcare facilities, aged care, and the community.

Since being established in 2011 we've put customers first. But our origins date back to 1910 in Chicago when Medline manufactured butcher's aprons, and soon after, surgical gowns. And the rest is history. So, with Medline you get a local team that understand you, backed by a global healthcare company that's been doing it for generations.

Together, advancing healthcare.

Table of contents

Introduction	3
Sterilisation wrap	4
Sterilisation accessories.....	7
Sterilisation pouches and rolls	8
Gloves.....	9





Introduction

Your CSSD revolves around your OR. We get that - and that's why Medline's range of products are designed to assist you through each step of the process: from packaging to storage, ensuring you have the right tools to aid in sterilisation and proper safekeeping has never been more important.

Packaging is the critical bridge between preparation and sterilisation. That's why from wraps, to trays and accessories, our easy-to-use products are designed to simplify your work, increase efficiency and help you assemble sterilised instrument trays with confidence.

Sterilisation wrap

Choose the right Gemini wrap for you - whether that be two single sheets wrapped separately, or bonded sheets to wrap simultaneously.

Which ever you choose, both will allow sterilization of the enclosed medical devices, whilst also maintaining sterility until use.



100% polypropylene construction provides high-level barrier protection against fluid and particulates



Strong for resistance against holes and tears, with softer, smooth finish easier to handle and fold



Drapeable with less memory, supporting operating theatre staff when opening



Versatility, indicated for use in all major sterilisation cycles:

- Pre-vacuum Steam - STERRAD®
- Gravity Steam - Steris V-Pro®
- Ethylene Oxide (ETO)

Bonded wrap: The right Gemini for the job

- Ultrasonically fused at the seam, to provide exceptional strength and improved efficiency
- Comprised of 5 layer SSMMS polypropylene
- Dual blue/pink colouring to differentiate specific instrument sets
- Facilitates easier detection of perforations and strike through



Gemini Bonded Wrap is available in a wide range of sizes and thickness for all types of instrument trays.



What are you sterilising today?**We have a wrap for that.**

From light weight to super weight, we've got more ways for you to wrap with confidence.


SINGLE PLY

	Light weight	Regular weight	Medium weight	Heavy weight	Super heavy weight
Recommended maximum weight	2.7kg 47gsm	4kg 54gsm	5.9kg 61gsm	7.7kg 71gsm	11.3kg 75gsm
18" x 18" (45.7 x 45.7 cm)	GEM1118-EU (1000/Cs)	GEM2118-EU (1000/Cs)	-	-	-
20" x 20" (50.8 x 50.8 cm)	-	GEM2120-EU (500/Cs)	-	-	-
24" x 24" (61 x 61 cm)	GEM1124-EU (500/Cs)	GEM2124-EU (500/Cs)	GEM3124-EU (250/Cs)	GEM4124-EU (250/Cs)	-
30" x 30" (76.2 x 76.2 cm)	GEM1130-EU (300/Cs)	GEM2130-EU (300/Cs)	GEM3130-EU (200/Cs)	GEM4130-EU (250/Cs)	-
36" x 36" (91.4 x 91.4 cm)	GEM1136-EU (300/Cs)	GEM2136-EU (300/Cs)	GEM3136-EU (150/Cs)	GEM4136-EU (150/Cs)	GEM5136-EU (144/Cs)
40" x 40" (101.6 x 101.6 cm)	-	GEM2140-EU (250/Cs)	GEM3140-EU (150/Cs)	-	GEM5140-EU (96/Cs)
45" x 45" (114.3 x 114.3 cm)	-	-	GEM3145-EU (100/Cs)	GEM4145-EU (100/Cs)	-
48" x 48" (121.9 x 121.9 cm)	GEM1148-EU (250/Cs)	GEM2148-EU (100/Cs)	GEM3148-EU (50/Cs)	GEM4148-EU (50/Cs)	GEM5148-EU (48/Cs)
54" x 54" (137.2 x 137.2 cm)	-	GEM2154-EU (100/Cs)	GEM3154-EU (50/Cs)	-	-
54" x 72" (137.2 x 182.9 cm)	-	-	GEM3172-EU (50/Cs)	GEM4172-EU (50/Cs)	-

DOUBLE PLY

	Light weight	Regular weight	Medium weight	Heavy weight	Super heavy weight
Recommended maximum weight	2.7kg 47gsm	4kg 54gsm	5.9kg 61gsm	7.7kg 71gsm	11.3kg 75gsm
18" x 18" (45.7 x 45.7 cm)	GEM1118T-EU (500/Cs)	-	-	-	-
20" x 20" (50.8 x 50.8 cm)	GEM1120T-EU (500/Cs)	-	-	-	-
24" x 24" (61 x 61 cm)	GEM1124T-EU (250/Cs)	GEM2124T-EU (250/Cs)	GEM3124T-EU (100/Cs)	-	-
30" x 30" (76.2 x 76.2 cm)	GEM1130T-EU (150/Cs)	GEM2130T-EU (150/Cs)	-	GEM4130T-EU (100/Cs)	-
36" x 36" (91.4 x 91.4 cm)	GEM1136T-EU (150/Cs)	GEM2136T-EU (150/Cs)	GEM3136T-EU (75/Cs)	GEM4136T-EU (75/Cs)	-
40" x 40" (101.6 x 101.6 cm)	GEM1140T-EU (100/Cs)	GEM2140T-EU (100/Cs)	GEM3140T-EU (75/Cs)	-	GEM5140T-EU (48/Cs)
45" x 45" (114.3 x 114.3 cm)	-	-	GEM3145T-EU (50/Cs)	GEM4145T-EU (50/Cs)	GEM5145T-EU (48/Cs)
48" x 48" (121.9 x 121.9 cm)	GEM1148T-EU (100/Cs)	GEM2148T-EU (50/Cs)	GEM3148T-EU (30/Cs)	GEM4148T-EU (30/Cs)	GEM5148T-EU (24/Cs)
54" x 54" (137.2 x 137.2 cm)	-	-	GEM3154T-EU (30/Cs)	GEM4154T-EU (30/Cs)	GEM5154T-EU (24/Cs)
54" x 72" (137.2 x 182.9 cm)	-	-	-	GEM4172T-EU (30/Cs)	GEM5172T-EU (24/Cs)

For more information refer to the
sterilisation wrap brochure



Sterilisation tray liner

- Can be used in either wrapped instrument trays or within containers
- Assists in wicking away any excess moisture remaining in post sterilisation

Item number	Dimensions	Pkg
WGN400	51 x 64cm	450/Cs



WGN400

Sterilisation transport tray

- Rigid plastic with perforations and ridged bottom to permit air flow
- Used to store and transfer wrapped instrument sets

Item number	Dimensions	Pkg
GEMTRAY3CE	30 x 56cm	10/box




GEMTRAY3CE

Disposable tray corner protector

- Foam material
- Placed on the corners of the sterilisation tray
- Single use

Item number	Dimensions	Pkg
GEMCORNER	15.2 x 15.2cm	100/Cs



GEMCORNER



Sterilisation pouches

- Sterilisation pouch for sharps and other instruments
- For use in steam or EO sterilisation
- Pre-sealed on three edges, self-seal last edge
- Colour change indicators for steam or EO sterilisation identification

Item number	Dimensions	Pkg
MPP100520GS	8.9 x 22.9cm	200ea/box, 2,000ea/Cs
MPP100525GS	13.3 x 25.4cm	200ea/box, 2,000ea/Cs
MPP100545GS	19 x 33cm	200ea/box, 1,000ea/Cs
MPP100546GS	20.3 x 40.6cm	200ea/box, 1,000ea/Cs
MPP100550GS	12.7 x 38.1cm	200ea/box, 1,000ea/Cs
MPP100555GS	30.5 x 38.1cm	200ea/box, 1,000ea/Cs
MPP100558GS	30.5 x 45.7cm	200ea/box, 1,000ea/Cs
MPP100565GS	5.9 x 55.9cm	200ea/box, 1,000ea/Cs

Tyvek pouches

Item number	Dimensions	Pkg
MPP100650N (Tyvek)	8.9 x 22.9cm	200ea/box, 2,000ea/Cs
MPP100655N (Tyvek)	13.3 x 26cm	200ea/box, 2,000ea/Cs

Sterilisation rolls

- Sterilisation roll for use in steam auto claves and EO
- Can be cut to different lengths for versatility
- Ideal for odd-sized, long or bulky objects
- Requires heat sealing on both ends
- Maintains sterility up to 1 year post sterilisation
- Colour change indicators for steam or EO sterilisation identification



Item number	Dimensions	Pkg
MPP200410R	7.6cm x 200m	1/Cs
MPP200410S	7.6cm x 30.5m	8/Cs
MPP200415R	10.1cm x 200m	1/Cs
MPP200415S	10.1cm x 30.5m	16/Cs
MPP200420R	15.2cm x 200m	1/Cs
MPP200420S	15.2cm x 30.5m	8/Cs
MPP200422R	20.3cm x 200m	1/Cs
MPP200422S	20.3cm x 30.5m	4/Cs
MPP200425R	25.4cm x 200m	1/Cs
MPP200425S	25.4cm x 30.5m	4/Cs

Gloves

Our SensiCare Pro 400 gloves are ideal for the protection of professionals who use harsh chemicals.


- Range of sizes from Small to X-Large
- Colour coded packaging for easier size detection
- Uniform box and case measurements across all sizes
- Flat-layer packing to prevent overlapping within the box extended cuff
- Textured grip to enhance tactile sensitivity in both wet and dry situations

Item number	Size	Pkg
PRO400S	Small	50ea/box, 10box/Cs
PRO400M	Medium	
PRO400L	Large	
PRO400XL	X-Large	



Size	Small	Medium	Large	X-large
	Orange	Blue	Green	Black





For more information, contact your Medline Account Manager or call 1800 110 511.

**Together,
advancing
healthcare**



Medline International Two Australia Pty Ltd
2 Fairview Place, Marsden Park NSW 2765
1800 110 511
medline.com.au