

Gemini Sterilisation Wraps



**Versatile
barrier
protection**

About Medline


We're a global healthcare company dedicated to being easy to partner with, solving problems, and providing good value. In Australia, from our network of warehouses nationally, we manufacture and distribute an expansive portfolio of products for hospitals, healthcare facilities, aged care, and the community.

Since being established in 2011 we've put customers first. But our origins date back to 1910 in Chicago when Medline manufactured butcher's aprons, and soon after, surgical gowns. And the rest is history. So, with Medline you get a local team that understand you, backed by a global healthcare company that's been doing it for generations.

Together, advancing healthcare.

Table of contents

Single ply.....	3
Double ply.....	4



Single ply

With Gemini Wrap, you have a choice: fold two wraps at the same time, or sequentially.

- **Made with SMS fabric:** Gemini offers protection against fluid and particles
- **Low linting:** Gemini wrap is naturally low linting to help maintain the integrity of the sterile trays
- **Drapeable:** Gemini wrap is very soft and easy to fold, making it simple to work with
- **Adhesive properties:** The surface is compatible with sterilisation tape
- **Strong:** Gemini wrap is resistant to tears and punctures



	Light weight	Regular weight	Medium weight	Heavy weight	Super weight
Recommended maximum weight	2.7kg 47gsm	4kg 54gsm	5.9kg 61gsm	7.7kg 71gsm	11.3kg 75gsm
18" x 18" (45.7 x 45.7 cm)	GEM1118-EU (1000/cs)	GEM2118-EU (1000/cs)	-	-	-
20" x 20" (50.8 x 50.8 cm)	-	GEM2120-EU (500/cs)	-	-	-
24" x 24" (61 x 61 cm)	GEM1124-EU (500/cs)	GEM2124-EU (500/cs)	GEM3124-EU (250/cs)	GEM4124-EU (250/cs)	-
30" x 30" (76.2 x 76.2 cm)	GEM1130-EU (300/cs)	GEM2130-EU (300/cs)	GEM3130-EU (200/cs)	GEM4130-EU (250/cs)	-
36" x 36" (91.4 x 91.4 cm)	GEM1136-EU (300/cs)	GEM2136-EU (300/cs)	GEM3136-EU (150/cs)	GEM4136-EU (150/cs)	GEM5136-EU (144/cs)
40" x 40" (101.6 x 101.6 cm)	-	GEM2140-EU (250/cs)	GEM3140-EU (150/cs)	-	GEM5140-EU (96/cs)
45" x 45" (114.3 x 114.3 cm)	-	-	GEM3145-EU (100/cs)	GEM4145-EU (100/cs)	-
48" x 48" (121.9 x 121.9 cm)	GEM1148-EU (250/cs)	GEM2148-EU (100/cs)	GEM3148-EU (50/cs)	GEM4148-EU (50/cs)	GEM5148-EU (48/cs)
54" x 54" (137.2 x 137.2 cm)	-	GEM2154-EU (100/cs)	GEM3154-EU (50/cs)	-	-
54" x 72" (137.2 x 182.9 cm)	-	-	GEM3172-EU (50/cs)	GEM4172-EU (50/cs)	-

Double ply

Gemini bonded wrap can provide you greater confidence in your sterilisation process. Two sheets electronically sealed together provide exceptional strength and improve efficiency.

Barrier protection: Strong dual-layer, 100% polypropylene construction provides barrier protection against fluid and particles

Easy to use: Drapable with less memory when opening. Soft, smooth finish is easy to handle and fold

Versatile: Medline's Gemini Bonded Wrap is for use in all major sterilisation cycles, including:


- Pre-Vacuum Steam Cycles
- Gravity Steam Cycles
- Ethylene Oxide (ETO) Sterilisation
- STERRAD® Sterilisation
- STERIS V-PRO® Sterilis

Our dual colour bonded blue/pink allows for:

- Quicker identification of loan sets
- Easier detection of tears, punctures and strike through



	Light weight	Regular weight	Medium weight	Heavy weight	Super weight
Recommended maximum weight	2.7kg 47gsm	4kg 54gsm	5.9kg 61gsm	7.7kg 71gsm	11.3kg 75gsm
18" x 18" (45.7 x 45.7 cm)	GEM1118T-EU (500/cs)	-	-	-	-
20" x 20" (50.8 x 50.8 cm)	GEM1120T-EU (500/cs)	-	-	-	-
24" x 24" (61 x 61 cm)	GEM1124T-EU (250/cs)	GEM2124T-EU (250/cs)	GEM3124T-EU (100/cs)	-	-
30" x 30" (76.2 x 76.2 cm)	GEM1130T-EU (150/cs)	GEM2130T-EU (150/cs)	-	GEM4130T-EU (100/cs)	-
36" x 36" (91.4 x 91.4 cm)	GEM1136T-EU (150/cs)	GEM2136T-EU (150/cs)	GEM3136T-EU (75/cs)	GEM4136T-EU (75/cs)	-
40" x 40" (101.6 x 101.6 cm)	GEM1140T-EU (100/cs)	GEM2140T-EU (100/cs)	GEM3140T-EU (75/cs)	-	GEM5140T-EU (48/cs)
45" x 45" (114.3 x 114.3 cm)	-	-	GEM3145T-EU (50/cs)	GEM4145T-EU (50/cs)	GEM5145T-EU (48/cs)
48" x 48" (121.9 x 121.9 cm)	GEM1148T-EU (100/cs)	GEM2148T-EU (50/cs)	GEM3148T-EU (30/cs)	GEM4148T-EU (30/cs)	GEM5148T-EU (24/cs)
54" x 54" (137.2 x 137.2 cm)	-	-	GEM3154T-EU (30/cs)	GEM4154T-EU (30/cs)	GEM5154T-EU (24/cs)
54" x 72" (137.2 x 182.9 cm)	-	-	-	GEM4172T-EU (30/cs)	GEM5172T-EU (24/cs)



For more information, contact your Medline Account Manager or call 1800 110 511.

**Together,
advancing
healthcare**



Medline International Two Australia Pty Ltd
2 Fairview Place, Marsden Park NSW 2765
1800 110 511
medline.com.au

We reserve the right to correct any errors that may occur within this brochure. © 2023 Medline Industries, LP. All rights reserved.
Medline is a registered trademark of Medline Industries, LP. OR015/Aug23/v4.
All products listed are class I non-sterile medical devices. These products are intended to be used by healthcare professionals.
Before use, consult instructions and precautions on the corresponding labelling.