

# OR Necessities

Separate sterile pack components



**Standardisation  
and efficiency**

# About Medline

We're a global healthcare company dedicated to being easy to partner with, solving problems, and providing good value. In Australia, from our network of warehouses nationally, we manufacture and distribute an expansive portfolio of products for hospitals, healthcare facilities, aged care, and the community.

In Australia, since being established in 2011 we've put customers first. But our origins date back to 1910 in Chicago when Medline manufactured butcher's aprons, and soon after, surgical gowns. And the rest is history. So, with Medline you get a local team that understand you, backed by a global healthcare company that's been doing it for generations.

**Together, advancing healthcare.**

## Table of contents

Markers and rulers .....	4
Needle counters.....	5
Light handle covers .....	6
Accessories.....	6
Procedure prep.....	7
Fluid waste bowl.....	7



# Introduction

## Achieve consistency in off-the-shelf components

Medline's OR Necessities portfolio includes the sterile components found in your surgical packs. Separate pack components will help you to:



### Improve standardisation

- Easy product access – gain easy access to replacements for components contaminated during procedures or when additional components are required.
- Fewer SKUs – we review items on your shelf and in your packs, recommending ways to improve consistency and eliminate waste.



### Minimise risk

- Consistent practices across pack components and sterile shelf stock.
- Supporting sharps safety by meeting Australian healthcare standards for sharps disposal.



### Contain costs

- Stock only necessary pack components and reduce waste related costs.
- Increase set-up efficiencies and reduce associated labour costs.

Medline strives to be a true partner with our customers and clinicians, helping to find new ways to help our customers improve and optimise their safety practices. Medline's OR consumables product line is designed to make best-practice techniques more intuitive for clinicians. You can add these components in your packs and have them separately on your shelf. They're easy to implement and require little procedural change.



# Intraoperative patient prep

Our full line of skin markers, needle counters, instrument mats and other accessories can minimise risk and support patient and clinical safety when used within guidelines and the surgical safety checklist.

## Markers and rulers

### Skin markers

Medline's skin markers contain medical-grade, non-irritating gentian violet ink.

- Large ink reservoir for longer product life
- Constant ink flow for clear marking
- Choice of regular and fine-tip styles
- Available with ruler and labels

Item number	Description	Pkg	Non-sterile component
DYNJSM01	Sterile, regular tip with ruler	50ea/cs	MD2001BNR
DYNJSM02	Sterile, regular tip with ruler and 6 labels	50ea/cs	MD2001BNRL
DYNJSM03	Sterile, fine tip with ruler	50/cs	



### Utility marker

- Black, non-smearing ink
- Can be used on most surfaces for labelling

Item number	Description	Pkg
DYNJSM06	Sterile, utility tip with ruler	50ea/cs





## Needle counters

- Sturdy plastic container keeps track of used disposable surgical sharps during procedure
- Available in foam block, adhesive, foam strip and magnet configurations for clinician preference
- Adhesive tabs secure the counters to workspaces like mayo table covers
- Counts range from 10 to 100

Item number	Description	Pkg	Non-sterile component
NCF10MCE	10 count, foam block with magnet	16ea/cs	
DYNJNC15FSMA	15 count, foam strip with magnet	168ea/cs	DYNJNC15FSMABNS
GSNCF20MBRCE	20 count, foam block with magnet	25/cs	49814AU
NC20MMCEA	20 count, double magnet	168/cs	NC10201MM

Item number	Description	Pkg	Non-sterile component
NC30MMBRCE	30 count, double magnet	25ea/cs	
NC30FSMBRA	30 count, foam strip with magnet	50ea/cs	NC3060FMGY
NC40FBRA	40 count, foam block with magnet	50/cs	49255AU
DYNJNC100FA	100 count, double magnet with foam block	64/cs	
NDLNESTCE	Needle nest	25/cs	35862



## Light handle covers

- Flexible, semi-rigid light handle covers are pushed onto the light handle
- Packed in a peel pouch either single or double pack
- Universal light handle covers fit on most light handles without using an adaptor



Item number	Description	Pkg
DYNJLHS1	Green light handle cover, single pack	80ea/cs
DYNJLHS2	Green, light handle cover, double pack	160ea/cs
16558A	Universal light handle cover	200/cs

## Accessories

### Magnetic instrument mat

- The magnetic instrument mat allows for the quick access and protection of surgical instruments.
- Flexible design conforms to the contours of the patient's body
- Non-slip foam backing helps prevent shifting over of surgical drape
- Magnets reduce the risk of instruments falling off the sterile field



### Anti-skid instrument mat

- Brightly coloured mat provides a very visible neutral zone for hands-free passing of sharps
- Textured surface grips instruments and helps prevent them from falling off the sterile field
- Can be placed on the sterile field or on the mayo stand
- Has two adhesive tapes on the bottom



Item number	Description	Pkg	Non-sterile component
MDL1CE	41 x 51 cm	20ea/cs	MDL1CECS
GSASM1CE	25 x 41 cm	80ea/cs	82658

## Procedure prep

### Anti-fog solution

- Solution helps surgeons see clearly, minimising fog forming on cameras and scopes
- Adhesive-backed sponge included for easy application of the solution



### Scrub brush

- Medline's scrub brushes help to achieve maximum cleaning with minimal skin irritation
- The flexible design fits comfortably in most hands
- Nail cleaner included



Item number	Description	Pkg	Non-sterile component
NONFB100	Anti-fog solution with sponge	20ea/cs	NONFB100CS
DYND371603	Dry scrub brush	80ea/cs	
DYNDBRKT	Scrub brush bracket	1ea	



## Fluid waste bowl

- Medline's fluid waste bowl assists with efficient disposal of fluids, improving patient and clinician safety by minimising the risk of splash back

Item number	Description	Pkg	Non-sterile component
GSSPLASH1CE	Fluid waste bowl 500ml, gold	24ea/cs	95666







**Medline offers a full range of sterile custom procedure packs - with over 2 million\* produced each year from our Marsden Park facility in Sydney's west.**

\* Figure calculated on 2022 Medline custom procedure pack sales.

**For more information, contact your Medline Account Manager or call 1800 110 511.**

**Together,  
advancing  
healthcare**



**Medline International Two Australia Pty Ltd**  
2 Fairview Place, Marsden Park NSW 2765  
1800 110 511  
[medline.com.au](https://www.medline.com.au)

Skin markers, light handle covers and instrument mats are class I sterile medical devices intended to be used by healthcare professionals.  
Anti-fog solution NONFB100 is a class IIa sterile medical device intended to be used by healthcare professionals.  
Before use, consult instructions and precautions on the corresponding labelling.  
We reserve the right to correct any errors that may occur within this brochure. © 2022 Medline Industries, LP. All rights reserved. OR057/Jan23/v1.