

QuickSuite® OR table kits and sheets



The complete single use solution for your
operating theatre



**Improve
efficiency
Promote
protection**

About Medline


We're a global healthcare company dedicated to being easy to partner with, solving problems, and providing good value. In Australia, from our network of warehouses nationally, we manufacture and distribute an expansive portfolio of products for hospitals, healthcare facilities, aged care, and the community.

Since being established in 2011 we've put customers first. But our origins date back to 1910 in Chicago when Medline manufactured butcher's aprons, and soon after, surgical gowns. And the rest is history. So, with Medline you get a local team that understand you, backed by a global healthcare company that's been doing it for generations.

Together, advancing healthcare.

Table of contents

QuickSuite®	3
Sahara® absorbent sheets.....	4
Sahara® table sheets.....	5
QuickSuite® Stock Kits.....	6
Ordering information & QuickSuite® accessories....	7
Liqui-Loc™	8



QUICKSUITE®

Helps your OR operate smarter

Medline's QuickSuite® range includes disposable operating table sheets and kits to elevate OR set up and turnover. QuickSuite® aims to improve some of your hospital and theatre challenges by helping you to:



Improve efficiency

- Reduce the time it takes to turn over your OR with our complete kits.
- Less touchpoints for staff during set up and clean up allows for more time in theatre.



Reduce risk

- Lower the risk of cross contamination with Sahara® absorbent OR table sheets.
- All our QuickSuite® sheets are purposely designed to be single use, making it easier to dispose of contaminants.



Lower costs


- Eliminate unnecessary waste with stock kits to avoid associated costs of reusable linens.

Sahara® absorbent table sheets

Switch to disposable to protect skin health

Help reduce the risk of moisture assisted skin damage with Sahara® absorbent table sheets

- 1 Quilted for smooth surface and patient comfort
- 2 Core polymer technology wicks moisture away maintaining skin integrity
- 3 Polyethylene film layer for maximum barrier protection



- Over 3.7L/m² of fluid absorbency
- Available in two sizes 102cm x 229cm (long) and 102cm x 152cm (short)
- Offered in kits which include draw/lift sheet, armboard covers and head covers
- Ideal for: long, high fluid procedures

Ordering information

Item no.	Description	Pkg
DYND4060SBE	Sahara® absorbent sheet; 102 x 152cm (short)	20/box
DYND4090SBE	Sahara® absorbent sheet; 102 x 229cm (long)	20/box

Sahara® table sheets

For your short simple procedures, our Sahara® economy table sheet and Sahara® protective sheet provides you with options for varying fluid output, rather than using one sheet for all procedures.

Sahara® economy OR table sheet

- One size for OR table use: 102 x 229cm
- Eliminates unnecessary linen laundry costs
- Ideal for: medium fluid procedures

Ordering information

Item no.	Description	Pkg
DYND4090E	Sahara® economy sheet 102 x 229cm	20/case



Sahara® protective OR table sheet

- Designed to protect your OR table and extend the life of your capital equipment
- Ideal for: short procedures

Ordering information

Item no.	Description	Pkg
DYND4090NA	Sahara® protective sheet 102 x 229cm	40/case 10/bag



QuickSuite® Stock Kits

Medline's QuickSuite® operating table stock kits contain an assortment of disposable components to be used prior a procedure to set up the operating theatre table, and afterwards to rapidly clear it for the next patient.

DYKQSUITEL2CEA



DYKQSUITEL3CEA



DYKQSUITEL4CEA




Item no.	Description	Sahara® Table Sheet 102 x 230cm	Draw/ Lift Sheet 102 x 152cm	Armboard Covers 72.5 x 33.5cm	Head Covers 38 x 39cm	Pkg
DYKQSUITEL2CEA	Sahara® absorbent sheet kit	•	•	-	-	24/box
DYKQSUITEL3CEA	Sahara® absorbent sheet kit	•	•	•	-	16/box
DYKQSUITEL4CEA	Sahara® absorbent sheet kit	•	•	•	•	16/box

Ordering information

Item no.	Description	Dimensions	Draw/Lift Sheet	Armboard Covers	Head Covers	Pkg
DYND4060SBE	Sahara® absorbent sheet (short)	102 x 152cm	-	-	-	20/box
DYND4090SBE	Sahara® absorbent sheet (long)	102 x 229cm	-	-	-	20/box
DYND4090E	Sahara® economy sheet	102 x 229cm	-	-	-	20/box
DYND4090NA	Sahara® protective sheet	102 x 229cm	-	-	-	40/box
DYKQSUITEL2CEA	QuickSuite® stock kit	102 x 229cm (DYND4090SBE)	•	-	-	24/box
DYKQSUITEL3CEA	QuickSuite® stock kit	102 x 229cm (DYND4090SBE)	•	•	-	16/box
DYKQSUITEL4CEA	QuickSuite® stock kit	102 x 229cm (DYND4090SBE)	•	•	•	16/box

QuickSuite® accessories



NON24398
Arm board covers with strap

Item no.	Description	Dimensions	Pkg
NON24398	Armboard cover with strap	33 x 72cm	100/cs
QSACOVER	Armboard cover strap	33 x 72cm	100/cs

*supplied in pairs



ORS01
Body straps

Item no.	Description	Dimensions	Pkg
ORS01	Disposable body strap	10 x 152cm	24/cs

*supplied in pairs

Liqui-Loc™

Also looking to protect your OR floor?

Medline's Liqui-Loc™ floor mats come in a range of sizes and absorbent capabilities to help prevent slips, trips and falls.


Ideal for:


- Reducing spill-related risks in operating theatres, exam room, scrub sinks
- Procedures where large volumes of fluid could be generated (i.e. arthroscopy, gynaecology, urology)



- Durable, non linting top layer
- Cold-form cellulose and superabsorbent polyester fibres
- Waterproof, anti-slip coloured backing

Item no.	Description	Dimensions	Absorbency	Pkg
ORLLPBS	Liqui-Loc™, blue, small	71 x 101cm	2.2L	13/cs
ORLLPBM	Liqui-Loc™, blue, medium	71 x 142cm	3.0L	10/cs
ORLLPPS	Liqui-Loc™, pink, small	71 x 101cm	4.6L	11/cs
ORLLPPL	Liqui-Loc™, pink, large	71 183cm	8.4L	6/cs





For more information, contact your Medline Account Manager or call 1800 110 511.

**Together,
advancing
healthcare**



Medline International Two Australia Pty Ltd
2 Fairview Place, Marsden Park NSW 2765
1800 110 511
medline.com.au