

ACCUtouch®

Non-sterile vinyl - powder-free - synthetic examination gloves



Reference(s)

MDS192073 (XS)
MDS192074 (S)
MDS192075 (M)
MDS192076 (L)
MDS192077 (XL)

Specifications

Material: Polyvinyl Chloride

Cuff: Beaded

Latex free: Yes

Accelerator: None

Grip:



Colour:



Thickness

Fingertip: 0.07 mm

Palm: 0.07 mm

Cuff: 0.06 mm

Dimensions

Box: 240 mm W x 120 mm L x 60 mm H

Case: 256 mm W x 315 mm L x 254 mm H

Packaging

100 gloves per box

10 boxes per case



Size

X-small

Small

Medium

Large

X-large

Yellow

Orange

Blue

Green

Black

Technical sheet

ACCUtouch®

Primary material	Polyvinyl Chloride
Colour	Clear
Grip	Smooth
Cuff design	Beaded
Latex free	Yes
Accelerator	None
Thickness	Finger tip 0.07 mm Palm 0.07 mm Cuff 0.06 mm
Dimensions	Box: 240 mm W x 120 mm L x 60 mm H Case: 256 mm W x 315 mm L x 254 mm H
Packaging	100 gloves per box, 10 boxes per case
Physical properties	EN 455-1 Freedom From Holes: 1.5 AQL EN 455-2 Dimensions: Length: 250 mm Width: 93 mm Force at Break: Before challenge: 4.7 N After challenge: 3.8 N EN 455-3 Powder Content: 0.12 milligrams per glove EN 455-4 Shelf Life Properties: Meets requirements and testing for shelf-life determination of 35 months
Chemical resistance	ACCUtouch® Synthetic exam gloves resistance to common chemicals has been assessed in accordance to EN 16523-1. Results can be obtained upon request.
Viral penetration	Tested and passed, in accordance with F1671
Regulations and quality standards	Class 1 non-sterile Medical Device per Regulation (EU) 2017/745. Complies with standards EN 455-1, 2, 3 and 4. Medline manufacturing locations are certified to EN ISO 13485.
Storage recommendations	Protect from freezing. Avoid excessive heat. Keep dry. Product should be shielded from direct sunlight, fluorescent lighting, X-rays, moisture and ozone.



For more information, contact your Medline Account Manager or call 1800 110 511.

**Together,
advancing
healthcare**



Medline International Two Australia Pty Ltd
2 Fairview Place, Marsden Park NSW 2765
1800 110 511
medline.com.au